



MEDICATION DELIVERY VIA CADD-PRIZM

GATHER SUPPLIES

Medication
IV Tubing
Battery - 9 Volt
Coin

PREPARATION PROCEDURE

- Wash Hands
- When using cassette & medication bag:
 - Open Tubing & Cassette from package
 - Remove Tip from End
 - SPIKE & Insert Tubing End into medication bag

INSTALLING BATTERY

- Press the button on the door
- Match + and –
(if battery is in backward the display will remain blank with no sound)

ATTACHING CASSETTE

- Attach Hook Hinge First
- Place upright on Flat Surface
- Push down & turn coin counter clockwise

PRIMING TUBING

- Press **OPTIONS**
- Press **ENTER**
- Press **Y** and hold until primed
- Press **NEXT** when completed

Starting Pump/Medication Infusion

- Take Cap off the end of the Tubing
- Wash White Cap with alcohol
- Push-Twist onto your White Connector
- Press **Stop/Start** Button
- Press **Y** to answer
 - Green light will Blink
 - ***“RUNNING”*** will appear
 - Fluid number will change

Stopping Pump

- Press **Stop/Start** button
- Press **Y** to answer
 - Amber light will Blink
 - ***“STOPPED”*** will appear

Changing Cassette

- Close Clamps
- Disconnect from your catheter
(**White Connector** stays with you)
- Using the coin – Turn clockwise
- Discard Medication Cassette as instructed
- Attach new Medication Cassette

ReSetting Your Pump

- Press **NEXT** to Reservoir Volume Screen
- Press **ENTER**
- Press **Y** to answer when “asked to Reset?”
- Prime as Instructed
- Start as Instructed