



**Approximate Cylinder Durations**  
**Continuous Flow – Measured in Hours**

C-Cylinder (Purple Band) ~240 liters					D-Cylinder (Black Band) ~425 liters				
PSI>	2000	1500	1000	500	PSI>	2000	1500	1000	500
LPM	Full	¾ Full	½ Full	¼ Full	LPM	Full	¾ Full	½ Full	¼ Full
1	4	3	2	1	1	7	5	3.5	1.5
2	2	1.5	1	Change	2	3.5	2.5	1.5	Change
3	1.5	1	Change	Change	3	2	1.5	1	Change
4	1	Change	Change	Change	4	1.5	1	Change	Change

E-Cylinder (Yellow Band) ~684 liters					M60-Cylinder (Orange Band) ~1700 liters				
PSI>	2000	1500	1000	500	PSI>	2000	1500	1000	500
LPM	Full	¾ Full	½ Full	¼ Full	LPM	Full	¾ Full	½ Full	¼ Full
1	11	8	5.5	3	1	28	21	14	7
2	5	3.5	2.5	1	2	14	11	7	3.5
3	3.5	2.5	2	Change	3	10	7	5	2.5
4	3	2	1.5	Change	4	7	5	3.5	1.5
5	2	1.5	1	Change					
6	1.5	1	Change	Change					

**Pulse Dose (OCR) – Measured in Hours**  
 (Values based on a Full Cylinder)

\*\*Cylinder duration may vary, depending on breaths per minute\*\*

Your Oxygen Flow Rate (LPM)	M6-Cylinder (Red Band) ~164 liters	C-Cylinder (Purple Band) ~240 liters	D-Cylinder (Black Band) ~425 liters	E-Cylinder (Yellow Band) ~682 liters
.5	16	24	42	69
.75	11	16	28	45
1	8	12	21	34
1.5	5.5	8	14	23
2	4	6	10	17
2.5	3	5	8	14
3	2.5	4	7	11.5
3.5	2	3.5	6	9
4	2	3	5	8.5

**Cylinder Sizes & Band Color**

